


## Typ T einziehbares curly Thermoelement Lead (JIS)




XE-7035-001

**TYPE-T (JIS)**



XE-7035-001

**TYPE-T  
(JIS)**



XE-7035-001

**TYPE-T (JIS)**

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Erhältlich in Typ T

Farbe kodiert zu JIS

2 Meter lang bei vollständiger Verlängerung

Beendet mit einem eingebauten Miniaturstecker

Polyurethan außen

### **Spezifikationen**

## **Specifications**

<b>Product Code</b>	XE-7035-001
<b>Cable Length</b>	2M fully extended
<b>Cable Type</b>	Polyurethane outer
<b>Connector Type</b>	Thermocouple Plug
<b>Connector Size</b>	Miniature
<b>Colour</b>	Brown
<b>Standards Met</b>	Colour code to JIS
<b>Thermocouple Type</b>	T
<b>Termination Type</b>	Miniature Thermocouple Plug