


## 4-adriger PFA-isolierter, verzinnter Kupferschirm, PRT-Sensorkabel / -draht

**10M**




**4 CORE**  
XE-2734-010

**25M**



**4 CORE**  
XE-2734-025

**10M**



**4 CORE**  
XE-2734-010

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**4-adriger PFA-isolierter, verzinnter Kupferschirm, PRT-Sensorkabel / -draht**

Hohe Temperatur bis 260C, extrudierte PFA-Teflon-Konstruktion. Wird hufig fr Pt100-Thermometer verwendet, auch mit anderen Typen wie Pt130, Pt500 und Pt1000 geeignet

**Leistungsbeschreibung****Specifications**

<b>General Description</b>	PFA insulated 7/0.2mm tin plated copper conductors, twisted cores, tin plated copper braided screen with PFA insulated outer jacket
<b>Cable Length</b>	10, 25m
<b>Cable Type</b>	PFA Teflon® construction
<b>Core / Strands</b>	7/0.2mm
<b>Insulation Material</b>	PFA Insulated
<b>Max. Temperature</b>	+260°C
<b>Min. Temperature</b>	-75°C
<b>Number of Wires</b>	4
<b>Standards Met</b>	Cores colour coded Red & White in accordance with IEC-751, black outer jacket