



## 6-adriges PFA-isolierter, verzinnter Kupferschirm, PRT-Sensorkabel / Draht



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**6-adriges PFA-isolierter, verzinnter Kupferschirm, PRT-Sensorkabel / Draht**

Hohe Temperatur bis 260°C, extrudiertes PFA Teflon® Konstruktion. Wird häufig für Pt100?-Thermometer verwendet, auch für andere Typen wie Pt130, Pt500 & Pt1000? geeignet

**Leistungsbeschreibung**
**Specifications**

<b>General Description</b>	PFA insulated 7/0.2mm tin plated copper conductors, twisted cores with nylon cord, tin plated copper braided screen with PFA insulated outer jacket
<b>Cable Length</b>	10, 25m, 50m
<b>Cable Type</b>	PFA Teflon® construction
<b>Core / Strands</b>	7/0.2mm
<b>Insulation Material</b>	PFA Insulated
<b>Max. Temperature</b>	+260°C
<b>Min. Temperature</b>	-75°C
<b>Number of Wires</b>	6
<b>Standards Met</b>	Cores colour coded Red & White in accordance with IEC-751, black outer jacket