



6-adriges PFA-isolierter, verzinnter Kupferschirm, PRT-Sensorkabel / Draht



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

6-adriges PFA-isolierter, verzinnter Kupferschirm, PRT-Sensorkabel / Draht

Hohe Temperatur bis 260°C, extrudiertes PFA Teflon® Konstruktion. Wird häufig für Pt100?-Thermometer verwendet, auch für andere Typen wie Pt130, Pt500 & Pt1000? geeignet

Leistungsbeschreibung
Specifications

General Description	PFA insulated 7/0.2mm tin plated copper conductors, twisted cores with nylon cord, tin plated copper braided screen with PFA insulated outer jacket
Cable Length	10, 25m, 50m
Cable Type	PFA Teflon® construction
Core / Strands	7/0.2mm
Insulation Material	PFA Insulated
Max. Temperature	+260°C
Min. Temperature	-75°C
Number of Wires	6
Standards Met	Cores colour coded Red & White in accordance with IEC-751, black outer jacket