



4-adriges PTFE-isolierter, versilberter Kupferschirm, PRT-Sensorkabel / -draht



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

4-adriges PTFE-isolierter, versilberter Kupferschirm, PRT-Sensorkabel / -draht

Kann verwendet werden, um 2-, 3- oder 4-Draht-RTD-Konfigurationen zu erweitern. Wird häufig für Pt100?-Thermometer verwendet, auch für andere Typen wie Pt130, Pt500 & Pt1000? geeignet

Leistungsbeschreibung**Specifications**

General Description	4 core PTFE insulated 7/0.2mm silver plated copper conductors, twisted cores, kapton tape, silver plated copper braided screen with PTFE insulated outer jacket
Cable Length	10, 25m
Cable Type	wrapped PTFE tape construction
Core / Strands	7/0.2mm
Insulation Material	Overall black PTFE insulation
Max. Temperature	+260°C
Min. Temperature	-75°C
Number of Wires	4
Standards Met	Cores colour coded 2 x Red / 2 x White to IEC-751 + black outer jacket