

## Bajonettkappen - Edelstahl



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

## Specifications

|                            |   |
|----------------------------|---|
| <b>General Description</b> | Bayonet caps provide a method of quick fitting into suitable adapters located in the process; this technique is widely used in plastics machinery |
| <b>Material</b>            | Stainless steel   |
| <b>Dimensions</b>          | 11.5mm ID to 12.5mm ID  |
| <b>Sizes</b>               | to suit spring 3/16" or 6.0mm   |