

## Bajonettadapter Edelstahl



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

## Specifications

<b>General Description</b>	Bayonet Adaptors available in various lengths and thread sizes
<b>Material</b>	Stainless steel
<b>Dimensions</b>	1/8" BSPP
<b>Length</b>	48mm or 64mm