

Solide gebohrte Thermotaschen



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Edelstahl-Thermotaschen in verschiedenen Größen sind für den Einsatz mit industriellen Temperaturfühler konzipiert, um die Ummantelung vor korrosiven Medien zu schützen und den Austausch der Sonde zu erleichtern, ohne den Prozess zu stören.

Ab Lager sind Prozessanschlüsse standardmäßig 1/2 Zoll BSPP, eine 1/8 Zoll BSPP Kompressionsform für den Durchmesser der Sonde ist erforderlich, um Sensor an Thermopocket anzuschließen.

Specifications

General Description	From stock, process connections are 1/2 inch BSPP as standard, a 1/8 inch BSPP compression fitting to suit diameter of probe is required to connect sensor to thermopocket.
Material	Stainless steel
Dimensions	303SS, 12mm Stem
Type & Stem Diameter	303SS, 12mm Stem
Stem Length	50mm, 100mm, 150mm
Probe Diameter	50mm, 100mm, 150mm stem length