


Solide gebohrte Hochdruck-Thermopocken



- NOTES:
1. Solid Bar Stock Drilled Thermowell
 2. Material 316/316L St.Stl
 3. Pocket Rated 5000PSI
 4. Material Certificate Required for each Batch

Dimensions in mm unless specified otherwise



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications

General Description	Solid Drilled High Pressure Thermopockets
Material	316/316L Stainless steel
Pocket Rated	5000 PSI
Dimensions	316SS Tapered (12mm at tip)
Length	100mm, 150mm, 200mm stem length
Type & Stem Diameter	316SS Tapered (12mm at tip)
Stem Length	100mm, 150mm, 250mm
Standards Met	Material Certificate supplied for each order