




Freiliegendes Thermoelement mit blankem Draht (0,025 mm Leiter) - Typ K



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Ultrafeines Thermoelement mit geschweiter Spitze Typ K fr optimale Empfindlichkeit.

Leistungsbeschreibung

Specifications

Product Code	XE-3517-001
General Description	Very fast response to a fraction of a second, Ideal for minimal invasive temperature sensing.
Sensor Type	Exposed Junction
Thermocouple Type	K
Cable Length	200mm
Cable Type	Fine gauge bare wire
Core / Strands	0.001" (0.025mm)
Max. Temperature	315°C