




Thermo Attachment Pads aus PTFE-Glasgewebe



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Selbstklebende Pads aus PTFE-Glasgewebe zur Lokalisierung von Thermoelementen mit geschweifter Spitze

Nenntemperatur -73 C bis 260 C

Padgre 13x18mm

Bogen enthlt 20 Blcke

Specifications

Product Code	XE-3599-001
Series	ATTACHMENT PADS