


Typ K IEC rechtwinklig - Kupfer-Scheibenspitze



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Edelstahlmantel mit geformtem Griff und 2 Meter Spiralverlängerungskabel, das in einem Miniatur-Thermoelementstecker abgeschlossen ist.

Specifications

Product Code	XE-3817-001
General Description	Surface temperature spring loaded thermocouple with right angled probe and copper disc tip
Thermocouple Type	K
Calibration	IEC
Cable Length	2m
Cable Termination	Fitted Plug
Max. Temperature	+600°C
Dimensions	Probe O8mm(17mmL) Tip x 60mm Hot Leg and O6mm x 180mm Cold Leg with O4.5mm copper disc.