


Typ K IEC-Rechtwinklige Oberfläche



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Edelstahlmantel mit geformtem Griff und 2 Meter gewickeltem Verlängerungskabel, das in einem Miniatur-Thermoelementstecker endet.

Specifications

Product Code	XE-3819-001
General Description	Surface temperature fast response thermocouple
Sensor Type	Right angled Probe
Thermocouple Type	K
Calibration	IEC
Cable Length	2m
Cable Termination	Fitted Plug
Max. Temperature	+900°C
Dimensions	Probe O12mm Tip x 40mm Hot Leg and O6mm x 180mm Cold Leg. Coiled element diameter 7mm
Tip	Ceramic Tip, Coiled Element