




LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

IEC-Feindraht-Scheiben-Thermoelement, PTFE - Typ K

XE-3515-001

FINE WIRE DISC THERMOCOUPLE TYPE K



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Kupferscheiben-Thermoelement mit Hochtemperatur-Klebestreifen (13x18mm x 5) zur schnellen und einfachen Befestigung an der Oberfläche/dem Prozess.

Specifications

General Description	Fast response – ideal for surface measurement
Sensor Type	Copper disc, 6.35mm diameter x 0.3mm thick
Thermocouple Type	K
Cable Length	1m
Cable Type	PTFE insulated twisted wire
Core / Strands	1/0.2mm
Termination Type	Bare tails
Max. Temperature	+200°C
Min. Temperature	-0°C
Standards Met	IEC-584-3