


Glasfaserisoliertes ANSI-Thermoelement mit Sichtkontakt und Miniaturstecker - Typ K



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Mssen Sie Ihren Temperatursensor kalibrieren?

Labfacility ist auf ISO17025 (UKAS)-Kalibrierstandards rckfhrbar und bietet die Wahl zwischen einer 3- oder 5-Punkt-rckfhrbaren Kalibrierung mit unserer hauseigenen Kalibriereinrichtung.

Bitte beachten Sie, dass Sie bezüglich der Lieferinformationen kontaktiert werden, wenn Sie Ihrer Bestellung eine Kalibrierungsoption hinzufügen.

Leistungsbeschreibung

Specifications

General Description	Fast response, welded exposed junction
Sensor Type	Good temperature resistance (not suitable in fluids)
Thermocouple Type	K
Thermocouple Junction	Exposed
Cable Length	300mm, 1, 2, 5 or 10 metre
Cable Type	glassfibre insulated flat pair cable
Core / Strands	1/0.3mm conductors
Cable Termination	Miniature Plug
Max. Temperature	+350°C (short periods up to 400°C)
Min. Temperature	-60°C
Tolerance	ANSI MC96.1 Class 2
Standards Met	ANSI MC96.1