


## Glasfaserisoliertes ANSI-Thermoelement mit Sichtschutz und blanken Schwänzen - Typ K



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Mssen Sie Ihren Temperatursensor kalibrieren?**

Labfacility ist auf ISO17025 (UKAS)-Kalibrierstandards rckfhrbar und bietet die Wahl zwischen einer 3- oder 5-Punkt-rckfhrbaren Kalibrierung mit unserer hauseigenen Kalibriereinrichtung.

**Bitte beachten Sie, dass Sie bezglich der Lieferinformationen kontaktiert werden, wenn Sie Ihrer Bestellung eine Kalibrierungsoption hinzufgen.**

**Leistungsbeschreibung****Specifications**

<b>General Description</b>	Fast response, welded exposed junction
<b>Sensor Type</b>	Good temperature resistance (not suitable in fluids)
<b>Thermocouple Type</b>	K
<b>Thermocouple Junction</b>	Exposed
<b>Cable Length</b>	300mm, 1, 2, 5 & 10 metre
<b>Cable Type</b>	glassfibre insulated flat pair cable
<b>Core / Strands</b>	1/0.3mm
<b>Cable Termination</b>	Bare tails
<b>Max. Temperature</b>	+350°C (short periods up to 400°C)
<b>Min. Temperature</b>	-60°C
<b>Tolerance</b>	ANSI MC96.1 Class 2
<b>Standards Met</b>	ANSI MC96.1