




## Platten-Blindbuchsen



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Wird zum Ausblenden von unbenutzten Plattenausschnitten verwendet, erhältlich für Miniaturplatten mit (FF)-Plattensockeln.

## Specifications

<b>General Description</b>	Used for blanking off unused panel cut-outs
<b>Connector Type</b>	FF
<b>Colour</b>	Black
<b>Sizes</b>	Miniature