




Platten-Blindbuchsen



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Wird zum Ausblenden von unbenutzten Plattenausschnitten verwendet, erhältlich für Miniaturplatten mit (FF)-Plattensockeln.

Specifications

General Description	Used for blanking off unused panel cut-outs
Connector Type	FF
Colour	Black
Sizes	Miniature