



**IEC-Blatt-Thermoelement, glasfaserisoliertes  
Kabel mit bergeflecht aus Edelstahl, das in  
blanken Schwitzen endet, Miniatur- oder  
Standardstecker, Typ K**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Leistungsbeschreibung****Specifications**

<b>General Description</b>	Thermocouple mounted on a stainless steel shim
<b>Thermocouple Type</b>	K or J
<b>Thermocouple Junction</b>	Junction grounded at tip for fast response
<b>Probe Diameter</b>	4.5mm
<b>Probe Length</b>	150mm
<b>Cable Length</b>	1, 2, 3, 4 or 5m
<b>Cable Type</b>	Glassfibre insulated lead with stainless steel overbraid
<b>Core / Strands</b>	7/0.2mm
<b>Termination Type</b>	Bare tails, Miniature or Standard Plug (with cable clamps)
<b>Max. Temperature</b>	+350°C
<b>Min. Temperature</b>	-60°C
<b>Dimensions</b>	25x13mm shim (0.2mm thick)