

Industrial Style Mineralinsulated Thermoelement mit KNE Klemmkopf - Typen K,J



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Benötigen Sie Ihren Temperatursensor kalibriert?

Nach ISO17025 (UKAS) Kalibrierstandards rückverfolgbar, bietet Labfacility die Wahl zwischen einer 3- oder 5-Punkt-Rückverfolgbaren Kalibrierung mit unserer hauseigenen Kalibrieranlage, sie können beliebige Temperaturpunkte zwischen -10°C & +200°C auswählen.

Bitte beachten Sie, dass, wenn Sie eine Kalibrierungsoption zu Ihrer Bestellung hinzufügen, dies 2-3 Tage Vorlaufzeit zu Ihrer Bestellung hinzufügt.

Specifications

General Description	Highly flexible, sheath can be bent/formed to suit many applications and processes
Sensor Type	Mineral insulated Thermocouple fitted with KNE terminal head
Thermocouple Type	Type K and J available
Thermocouple Junction	Insulated hot junction
Probe Diameter	6.0mm diameter
Sheath Material	310 stainless steel sheath
Head Type	KNE, Cast aluminium finish, screw type lid with chain
Gland	M16 x 1.5mm Plated brass gland fitted
Max. Temperature	+1100°C
Min. Temperature	-40°C
IP Rating	IP67