


# Hygienisches Pt100-Widerstandsthermometer, 1,5" RJT-Fassung mit montiertem Transmitter - 3-Leiter, 6,0 mm x 75 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Hygienisches Pt100-Widerstandsthermometer, 1,5-Zoll-RJT-Anschluss

## Specifications

<b>Product Code</b>	XE-5525-001
<b>General Description</b>	Ideal for hygienic processes including dairy, brewery and pharmaceutical applications.
<b>Sensor Type</b>	Pt100, 3 wire single element class A
<b>Head Type</b>	KNEF
<b>Connection Type</b>	Pt100, 3 wire single element class A
<b>Transmitter</b>	SEM-203p transmitter, ranged -50°C to +150°C
<b>Support Dia (D2)</b>	8.0mm
<b>Support Length (L2)</b>	50mm