


Hygienisches Pt100-Widerstandsthermometer, 1,5" RJT-Fassung mit montiertem Transmitter - 3-Leiter, 6,0 mm x 75 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Hygienisches Pt100-Widerstandsthermometer, 1,5-Zoll-RJT-Anschluss

Specifications

Product Code	XE-5525-001
General Description	Ideal for hygienic processes including dairy, brewery and pharmaceutical applications.
Sensor Type	Pt100, 3 wire single element class A
Head Type	KNEF
Connection Type	Pt100, 3 wire single element class A
Transmitter	SEM-203p transmitter, ranged -50°C to +150°C
Support Dia (D2)	8.0mm
Support Length (L2)	50mm