



Hygienisches Pt100-Widerstandsthermometer, 1,5" RJT-Anschluss mit montiertem Transmitter - 3-Leiter, 6,0 mm x 125 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

RJT-Verschraubungen (Ring Joint Type) werden oft als Standardverschraubung bezeichnet. Der Kompromiss einer O-Ring-Dichtung, die einen kleinen Spalt hinterlässt, in dem sich die Laufbuchse und der männliche Teil der Verbindung berlapen.

Ideal für hygienische Prozesse, einschließlich Molkerei-, Brauerei- und Pharmaanwendungen.

Leistungsbeschreibung

Specifications

Product Code	XE-5526-001
General Description	Ideal for hygienic processes including dairy, brewery and pharmaceutical applications.
Sensor Type	Pt100, 3 wire single element class A
Head Type	KNEF
Connection Type	Pt100, 3 wire single element class A
Transmitter	SEM-203p transmitter, ranged -50°C to +150°C
Support Dia (D2)	8.0mm
Support Length (L2)	50mm