



# Hygienisches Pt100-Widerstandsthermometer, 1,5" RJT-Anschluss mit montiertem Transmitter - 3-Leiter, 6,0 mm x 125 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

RJT-Verschraubungen (Ring Joint Type) werden oft als Standardverschraubung bezeichnet. Der Kompromiss einer O-Ring-Dichtung, die einen kleinen Spalt hinterlässt, in dem sich die Laufbuchse und der männliche Teil der Verbindung berlapen.

Ideal für hygienische Prozesse, einschließlich Molkerei-, Brauerei- und Pharmaanwendungen.

### Leistungsbeschreibung

## Specifications

<b>Product Code</b>	XE-5526-001
<b>General Description</b>	Ideal for hygienic processes including dairy, brewery and pharmaceutical applications.
<b>Sensor Type</b>	Pt100, 3 wire single element class A
<b>Head Type</b>	KNEF
<b>Connection Type</b>	Pt100, 3 wire single element class A
<b>Transmitter</b>	SEM-203p transmitter, ranged -50°C to +150°C
<b>Support Dia (D2)</b>	8.0mm
<b>Support Length (L2)</b>	50mm