

Hygienisches Pt100 Widerstandsthermometer, 1,5" Tri-Clamp-Befestigung



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Ideal für hygienische Prozesse, einschließlich Molkerei-, Brauerei- und Pharmaanwendungen.

Specifications

Sensor Type	Pt100, 3 or 4 wire single element class A
Probe Diameter	6mm
Probe Length	75mm OR 125mm
Head Type	KNEF
Transmitter	Optional SEM-203p transmitter, ranged -50°C to +150°C
Number of Wires	3