




LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Hygienisches Pt100 Widerstandsthermometer, 1,5" Tri-Clamp Anschluss ohne Transmitter - 4-Leiter, 6,0mm x 75mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Ideal fr hygienische Prozesse, einschließlich Molkerei-, Brauerei- und Pharmaanwendungen.

Leistungsbeschreibung

Specifications

Product Code	XE-5527-001
Sensor Type	Pt100, 3 or 4 wire single element class A
Probe Diameter	6.0mm
Probe Length	75mm
Head Type	KNEF
Transmitter	None
Number of Wires	3