

**Hygienisches Pt100 Widerstandsthermometer,  
1,5" Tri-Clamp-Anschluss mit montiertem  
Transmitter - 3-Leiter, 6,0 mm x 75 mm**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Ideal fr hygienische Prozesse, einschließlich Molkerei-, Brauerei- und Pharmaanwendungen.

### Leistungsbeschreibung

## Specifications

<b>Product Code</b>	XE-5529-001
<b>Sensor Type</b>	Pt100, 3 wire single element class A
<b>Probe Diameter</b>	6.0mm
<b>Probe Length</b>	75mm
<b>Head Type</b>	KNEF
<b>Transmitter</b>	SEM-203p transmitter, ranged -50°C to +150°C
<b>Number of Wires</b>	3