


**Hygienisches Pt100-Widerstandsthermometer,
1,5" Tri-Clamp-Anschluss mit montiertem
Transmitter - 3-Leiter, 6,0 mm x 125 mm**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Ideal fr hygienische Prozesse, einschließlich Molkerei-, Brauerei- und Pharmaanwendungen.

Leistungsbeschreibung

Specifications

Product Code	XE-5530-001
Sensor Type	Pt100, 3 wire single element class A
Probe Diameter	6.0mm
Probe Length	125mm
Head Type	KNEF
Transmitter	SEM-203p transmitter, ranged -50°C to +150°C
Number of Wires	3