

Pt100 3-Draht-Klasse B Widerstandsthermometer, KNE-Kopf mit Sender



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Benötigen Sie Ihren Temperatursensor kalibriert?

Nach ISO17025 (UKAS) Kalibrierstandards rückverfolgbar, bietet Labfacility die Wahl zwischen einer 3- oder 5-Punkt-Rückverfolgbaren Kalibrierung mit unserer hauseigenen Kalibrieranlage, sie können beliebige Temperaturpunkte zwischen -10°C & +200°C auswählen.

Bitte beachten Sie, dass, wenn Sie eine Kalibrierungsoption zu Ihrer Bestellung hinzufügen, dies 2-3 Tage Vorlaufzeit zu Ihrer Bestellung hinzufügt.

Specifications

General Description	Pt100 probe with KNE terminal head fitted with a transmitter
Sensor Type	Pt100 (100 Ohms @ 0°C) to IEC 751, Class B, 3 wire
Process Connection	1/2"BSP x 6mm
Calibration	IEC
Probe Diameter	6mm
Probe Length	150mm
Sheath Material	316 stainless steel
Head Type	KNE
Max. Temperature	+450°C (connection head @ 170°C)
Min. Temperature	-50°C
IP Rating	IP67
Accuracy	Class B
Transmitter	Yes
Range	Various ranges of transmitter available
Number of Wires	3