



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Pt100 3-Draht-Widerstandsthermometer der Klasse B, KNE-Kopf mit Messumformer - Messumformerbereich von 0-400 C



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Mssen Sie Ihren Temperatursensor kalibrieren?

Rckfhrbar auf ISO17025 (UKAS) Kalibrierstandards bietet Labfacility die Wahl zwischen einer 3- oder 5-Punkt-rckfhrbaren Kalibrierung mit unserer hauseigenen Kalibriereinrichtung, Sie knnen beliebige Temperaturpunkte zwischen -10 C und +200 C auswhlen.

Bitte beachten Sie, dass Sie bezglich der Lieferinformationen kontaktiert werden, wenn Sie Ihrer Bestellung eine Kalibrierungsoption hinzufgen.

Leistungsbeschreibung**Specifications**

Product Code	XE-5508-001
General Description	Pt100 probe with KNE terminal head fitted with a transmitter
Sensor Type	Pt100 (100 Ohms @ 0°C) to IEC 751, Class B, 3 wire
Process Connection	1/2"BSP x 6mm
Calibration	IEC
Probe Diameter	6.0mm
Probe Length	150mm
Sheath Material	316 Stainless Steel
Head Type	KNE
Max. Temperature	+450°C (connection head @ 170°C)
Min. Temperature	-50°C
IP Rating	IP67
Accuracy	Class B
Transmitter	0-400°C Range
Number of Wires	3