



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Pt100 4-Draht-Widerstandsthermometer der Klasse B, kompakter KNS-Kopf mit 1/2" BSP x 6 mm Anschluss mit 100 mm Sondenlänge



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Leistungsbeschreibung

Specifications

| | |
|---------------------------|--|
| Product Code | XE-5509-001 |
| Sensor Type | Pt100 (100 Ohms @ 0°C) to IEC 751, Class B, 4 wire |
| Process Connection | 1/2"BSP x 6mm |
| Probe Diameter | 6.0mm |
| Probe Length | 100mm |
| Sheath Material | 316 Stainless Steel |
| Head Type | KNS |
| Max. Temperature | +450°C (connection head @ 170°C) |
| Min. Temperature | -50°C |
| IP Rating | IP67 |
| Accuracy | Class B |
| Number of Wires | 4 |
| Series | Pt100 KNS |