


## PFA-Kabel aufgerumtes IEC-Thermoelement mit freiliegender Sperrschicht und montiertem Ministecker - Typen K,T



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Kabelaufgerumtes Thermoelement mit montiertem Miniatur-IEC-Stecker, PFA Twin Twisted oder Flat Pair.

Ein hochempfindliches, hochempfindliches Thermoelement mit feinem Draht und montiertem Miniaturstecker.

### Mssen Sie Ihren Temperatursensor kalibrieren?

Labfacility ist auf ISO17025 (UKAS) Kalibrierstandards rckfhrbar und bietet die Wahl zwischen einer rckfhrbaren 3- oder 5-Punkt-Kalibrierung mit unserer hauseigenen Kalibriereinrichtung, mit der Sie beliebige Temperaturpunkte zwischen -10 C und +200 C auswahlen knnen.

**Bitte beachten Sie, dass Sie, wenn Sie Ihrer Bestellung eine Kalibrierungsoption hinzufgen, bezglich der Lieferinformationen kontaktiert werden**

## Specifications

|                            |  |
|----------------------------|--|
| <b>General Description</b> | A complete, precision fine wire, high sensitivity thermocouple with a fitted miniature connector and cable tidy. |
| <b>Sensor Type</b>         | Exposed Junction   |
| <b>Cable Length</b>        | 1, 2m  |
| <b>Cable Type</b>          | PFA insulated wire, flat pair or twin twisted available  |
| <b>Termination Type</b>    | Miniature flat pin plug, colour coded in accordance IEC 584  |
| <b>Max. Temperature</b>    | Rated 220°C continuous   |
| <b>Tolerance</b>           | Class 1  |
| <b>Standards Met</b>       | IEC 584  |
| <b>Connector Type</b>      | Plug   |
| <b>Connector Size</b>      | Miniature  |