





# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## PFA-isoliertes ANSI-Thermoelement mit feiner Dicke und offener Sperrschicht (0,076 mm Leiter) - Typ K



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Mssen Sie Ihren Temperatursensor kalibrieren?**

Labfacility ist auf ISO17025 (UKAS)-Kalibrierstandards rckfhrbar und bietet die Wahl zwischen einer 3- oder 5-Punkt-rckfhrbaren Kalibrierung mit unserer hauseigenen Kalibriereinrichtung.

**Bitte beachten Sie, dass Sie bezglich der Lieferinformationen kontaktiert werden, wenn Sie Ihrer Bestellung eine Kalibrierungsoption hinzufügen.**

**Leistungsbeschreibung****Specifications**

<b>General Description</b>	Fine Gauge Exposed Welded Tip Thermocouple
<b>Sensor Type</b>	Fast response, welded tip exposed junction
<b>Thermocouple Type</b>	K
<b>Thermocouple Junction</b>	Exposed
<b>Cable Length</b>	300mm, 1m
<b>Cable Type</b>	Teflon® PFA insulated wires
<b>Core / Strands</b>	0.003" (0.076mm) conductors
<b>Cable Termination</b>	Bare tails
<b>Max. Temperature</b>	+260°C
<b>Min. Temperature</b>	-75°C
<b>Standards Met</b>	ANSI colour coded