



## 4 Kern Silikon Gummi isoliert, abgerundeter Abschnitt, PRT Sensorkabel / Draht



25M



**4 CORE**  
XE-2798-025

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**4 Kern Silikon Gummi isoliert, abgerundeter Abschnitt, PRT Sensorkabel / Draht**

Kann verwendet werden, um bis zu 4 Draht-RTD-Konfigurationen zu erweitern. Häufig verwendet mit Pt100"-Thermometern, aber geeignet für andere Typen Pt130/500/1000" usw.

**Spezifikationen****Specifications**

<b>General Description</b>	4 core silicone rubber insulated 7/0.2mm tin plated copper conductors, twisted together with silicone rubber outer jacket
<b>Cable Length</b>	10, 25m
<b>Cable Type</b>	Silicone rubber insulated
<b>Core / Strands</b>	7/0.2mm
<b>Insulation Material</b>	Silicone rubber insulated, tin plated copper conductors, twisted together with silicone rubber outer jacket
<b>Max. Temperature</b>	+200°C
<b>Min. Temperature</b>	-40°C
<b>Number of Wires</b>	4
<b>Standards Met</b>	Cores colour coded 2 x Red / 2 x White to IEC-751