




## Thermoelement-Kabel / -Draht (IEC) Typ E PVC-isoliertes Flachpaar



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Leistungsbeschreibung****Specifications**

|                            |   |
|----------------------------|---|
| <b>Cable Length</b>        | 10, 25, 50, 100, 200m                     |
| <b>Cable Type</b>          | Flat pair construction                    |
| <b>Core / Strands</b>      | Stranded conductors available, 7/0.2mm    |
| <b>Insulation Material</b> | PVC Insulated                             |
| <b>Max. Temperature</b>    | +105°C                                    |
| <b>Min. Temperature</b>    | -10°C                                     |
| <b>Standards Met</b>       | Colour code (cores & jacket) to IEC-584-3 |
| <b>Thermocouple Type</b>   | E   |
| <b>Tolerance</b>           | Class 2                                   |