

GM1311 Professionelles Hochleistungs-Digitales Lebensmittelthermometer



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

BITTE BEACHTEN SIE: Die Vorlaufzeit betrgt 7-10 Werktage bei der Auswahl der Kalibrierung

GM1311 Hochleistungs-Professionelles Digitales Lebensmittel-/Grill-/Fleischthermometer

Umgebungstemperatur: 0 bis 40 C (32 bis 104 F)

Luftfeuchtigkeit: 10-90% RH

Nennspannung: 3V

1 x Lebensmittel-Thermometer

1 x englisches Benutzerhandbuch

1 x CR2032

Leistungsbeschreibung

Specifications

Product Code	XE-4002-001
Accuracy	-10 - 100°C: +/-1°C / 100 - 200°C: +/-1.5°C
Resolution	0.1°C / °F
Colour	White
Dimensions	9.29 in x 1.34 in x 0.71 in (23.6 cm x 3.4 cm x 1.8 cm). Display dimension:44x18mm
Function #1	LCD display
Function #2	Temperature lock
Function #3	Max/Min measuring
Function #4	Conversion between Celsius degree and Fahrenheit degree
Function #5	Auto power off
Max. Temperature	300°C (572°F)
Min. Temperature	-50°C (-58°F)
Probe Length	149.5 x 6mm



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Weight	Weight: 0.95 oz (27 g)
Power Supply	CR2032 button battery x 1 (included)