

## GM1311 Professionelles Hochleistungs-Digitales Lebensmittelthermometer



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**BITTE BEACHTEN SIE: Die Vorlaufzeit betrgt 7-10 Werktage bei der Auswahl der Kalibrierung**

GM1311 Hochleistungs-Professionelles Digitales Lebensmittel-/Grill-/Fleischthermometer

Umgebungstemperatur: 0 bis 40 C (32 bis 104 F)

Luftfeuchtigkeit: 10-90% RH

Nennspannung: 3V

1 x Lebensmittel-Thermometer

1 x englisches Benutzerhandbuch

1 x CR2032

### Leistungsbeschreibung

## Specifications

<b>Product Code</b>	XE-4002-001
<b>Accuracy</b>	-10 - 100°C: +/-1°C / 100 - 200°C: +/-1.5°C
<b>Resolution</b>	0.1°C / °F
<b>Colour</b>	White
<b>Dimensions</b>	9.29 in x 1.34 in x 0.71 in (23.6 cm x 3.4 cm x 1.8 cm). Display dimension:44x18mm
<b>Function #1</b>	LCD display
<b>Function #2</b>	Temperature lock
<b>Function #3</b>	Max/Min measuring
<b>Function #4</b>	Conversion between Celsius degree and Fahrenheit degree
<b>Function #5</b>	Auto power off
<b>Max. Temperature</b>	300°C (572°F)
<b>Min. Temperature</b>	-50°C (-58°F)
<b>Probe Length</b>	149.5 x 6mm



# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

<b>Weight</b>	Weight: 0.95 oz (27 g)
<b>Power Supply</b>	CR2032 button battery x 1 (included)