


## IEC Velcro Thermoelement, PFA mit Miniaturstecker - Typ K,T



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Velcro-Thermoelemente können als schnelle und einfache Methode zur Messung der Temperatur um die Oberfläche eines Rohrs verwendet werden.

## Specifications

<b>Thermocouple Type</b>	K or T
<b>Cable Length</b>	1.5M (0.5M Inside Velcro)
<b>Cable Type</b>	PFA
<b>Core / Strands</b>	7/0.2mm
<b>Termination Type</b>	Miniature flat pin plug, colour coded in accordance IEC 584
<b>Max. Temperature</b>	+150°C
<b>Length</b>	Velcro: 500mm