


IEC-Zugmagnet-Thermoelement Typ K 3 kg, PFA-isoliertes Kabel mit Edelstahl-bergeflecht mit Miniatur- oder Standardstecker



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications

General Description	Rapid, reliable temperature sensing on ferrous surfaces
Thermocouple Type	K
Thermocouple Junction	Grounded measuring junction
Magnet	1" diameter (3kg pull)
Cable Length	1 or 2 metre
Cable Type	stainless steel over-braided PFA Teflon®
Cable Termination	Fitted Miniature or Standard Plug
Max. Temperature	+250°C (200°C Plug)