

**IEC Typ K 0,7 kg
Zugknopf-Magnet-Thermoelement, PTFE-isoliert
mit Miniatur- oder Standardstecker**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications

General Description	Ideal for ferrous metal surface measurement
Thermocouple Type	K
Thermocouple Junction	Grounded at tip
Magnet	12.5mm diameter button type magnet. Max pull: 0.7 KG
Cable Length	1 or 2 metres
Cable Type	PTFE insulated twin-twist wires
Cable Termination	Miniature or Standard Fitted Plug
Max. Temperature	+250°C (plug 200°C)