

## Mineralisiertes Dual-Thermoelement mit Standard-Duplex-Typ-K-Stecker



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Mssen Sie Ihren Temperatursensor kalibrieren?**

Labfacility ist auf ISO17025 (UKAS)-Kalibrierstandards rckfhrbar und bietet die Wahl zwischen einer 3- oder 5-Punkt-rckfhrbaren Kalibrierung mit unserer hauseigenen Kalibriereinrichtung.

**Bitte beachten Sie, dass Sie bezglich der Lieferinformationen kontaktiert werden, wenn Sie Ihrer Bestellung eine Kalibrierungsoption hinzufgen.**

**Leistungsbeschreibung****Specifications**

<b>General Description</b>	310 stainless steel sheath. Sheath can be bent/formed to suit many applications and processes.
<b>Sensor Type</b>	Mineral insulated Type 'K' with Duplex Thermocouple Plug
<b>Thermocouple Type</b>	K IEC
<b>Thermocouple Junction</b>	Insulated hot junction
<b>Termination Type</b>	Standard duplex plug termination (200°C)
<b>Max. Temperature</b>	1100°C
<b>Min. Temperature</b>	-40°C
<b>Standards Met</b>	Conforms to IEC 584 specification
<b>Connector Type</b>	Plug
<b>Connector Size</b>	Standard
<b>Length</b>	150, 250 or 500mm