

**Keramik Ofen Thermoelemente für hohe  
Temperatur (1400A, A°C) mit auminAfA¶sen  
PorzellanhAfA¼llen - Typen K,N**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Das 100 mm lange Trägerrohr aus Edelstahl befestigt den Keramikmantel am Industriekopf und ist in der Gesamtmantellänge enthalten.

## Specifications

<b>General Description</b>	Ceramic thermocouples for furnaces, kilns and other high temperature applications
<b>Sensor Type</b>	Made with rugged 20 gauge (0.8mm dia) thermocouple wire for long-term performance
<b>Thermocouple Type</b>	Types K and N with Aluminous Porcelain sheaths, rated to 1400°C
<b>Sheath Material</b>	12 mm diameter ceramic sheath
<b>Max. Temperature</b>	+1400°C
<b>IP Rating</b>	IP67 rated industrial protection head with cable entry and screw terminal block
<b>Support Tube</b>	21 mm O.D. support tube, 100 mm long