


Platin-Widerstand Pt100 Industriesonde mit Lagging Extension



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Benötigen Sie Ihren Temperatursensor kalibriert?

Nach ISO17025 (UKAS) Kalibrierstandards rückverfolgbar, bietet Labfacility die Wahl zwischen einer 3- oder 5-Punkt-Rückverfolgbaren Kalibrierung mit unserer hauseigenen Kalibrieranlage, sie können beliebige Temperaturpunkte zwischen -10°C & +200°C auswählen.

Bitte beachten Sie, dass, wenn Sie eine Kalibrierungsoption zu Ihrer Bestellung hinzufügen, dies 2-3 Tage Vorlaufzeit zu Ihrer Bestellung hinzufügt.

Specifications

Sensor Type	Pt100 (100 Ohms @ 0°C) to IEC 751
Insert	6.0mm diameter 316 stainless steel terminated with spring loaded ceramic connection block
Process Connection	1/2"BSP parallel
Cable Connection	M20 x 1.5mm cable entry (gland included)
Termination Type	IP67 aluminium alloy weatherproof KNE connection head
Max. Temperature	+450°C (connection head @ 170°C)
Min. Temperature	-50°C
Tolerance	Class B
Number of Wires	4
Probe Material	8.0mmO 316 stainless steel below 1/2"BSP parallel process connection 1/2" diameter x 75mm lagging length to underside of connection head