


GM520 Druckmanometer



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

[Klicken Sie hier.](#)

um die neueste Software herunterzuladen, die für den Betrieb dieses Produkts erforderlich ist.

Spezifikationen

Specifications

Product Code	XE-4011-001
Model	GM520 Pressure Manometer
Communications Protocol	USB communication
Accuracy	+/- 0.3% FSO (25°C)
Repeatability	+/- 0.2% (Maximum +/- 0.5% FSO)
Response	Typical 0.5 s
Linearity/hysteresis	+/- 0.29% FSO
Maximum Pressure	+/- 150kPa
Storage Temperature Range	-10 to 60°C
Ambient Operating Temperature	0 to 50°C
Measuring range	+/- 35kPa