

GM610 Holz FeuchtigkeitsmessgerÃt



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications

Product Code	XE-4013-001
Model	GM610 Wood Moisture Meter
Resolution	0.5%
Battery	2 x 1.5V AAA Battery
Humidity range	20%RH - 90%RH
Moisture Response Speed	<1s
Max. Temperature	+60°C
Min. Temperature	-10°C
Measuring range	2-40%/50%/60%/70%
Display	LCD
Series	Wood Moisture Meter