



## Fester Temperaturfühler fr Rennreifen



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

## Specifications

<b>Product Code</b>	XE-3891-001
<b>General Description</b>	Can be used with any digital thermometer with a type k mini input socket
<b>Sensor Type</b>	Fast response fixed racing tyre probe with a max depth of six millimetres
<b>Thermocouple Type</b>	K
<b>Calibration</b>	IEC
<b>Probe Length</b>	max depth of 6mm
<b>Cable Length</b>	1m
<b>Cable Termination</b>	Miniature Plug
<b>Max. Temperature</b>	+250°C
<b>Min. Temperature</b>	-50°C
<b>Features #1</b>	Junction Type: Grounded