

## 4-20mA Fernwandgehäuse, Pt100 Eingang mit PFA- oder SSOB-Kabel & PRT-Sensor



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Gehäuse in Allzweck-ABS geformt, bieten die externen Befestigungspunkte eine einfache und schnelle Möglichkeit, Sensoren direkt auf Kanäle zu montieren, ohne auf die Notwendigkeit von separaten Montageplatten.

### **Benötigen Sie Ihren Temperatursensor kalibriert?**

Nach ISO17025 (UKAS) Kalibrierstandards rückverfolgbar, bietet Labfacility die Wahl zwischen einer 3- oder 5-Punkt-Rückverfolgbaren Kalibrierung mit unserer hauseigenen Kalibrieranlage, sie können beliebige Temperaturpunkte zwischen -10°C & +200°C auswählen.

**Bitte beachten Sie, dass, wenn Sie eine Kalibrierungsoption zu Ihrer Bestellung hinzufügen, dies 2-3 Tage Vorlaufzeit zu Ihrer Bestellung hinzufügt.**

## **Specifications**

|                            |   |
|----------------------------|---|
| <b>General Description</b> | Compatible with Pt100 platinum resistance thermometers                |
| <b>Sensor Type</b>         | Surface Mounting for remote installation                              |
| <b>Cable Length</b>        | 1 or 2 metre  |
| <b>Cable Type</b>          | PFA screened or Fiberglass stainless steel overbraid available        |
| <b>Core / Strands</b>      | 7/0.2mm   |
| <b>Head Type</b>           | High Impact ABS Plastic Construction                                  |
| <b>Thread</b>              | Entry Threads M16   |
| <b>Transmitter</b>         | Yes   |
| <b>Range</b>               | -50°C to +150°C or 0°C to +400°C (either can be re-ranged if desired) |
| <b>Output</b>              | 4-20 mA output  |