

4-20mA Fernwandgehäuse, Pt100 Eingang mit PFA- oder SSOB-Kabel & PRT-Sensor



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Gehäuse in Allzweck-ABS geformt, bieten die externen Befestigungspunkte eine einfache und schnelle Möglichkeit, Sensoren direkt auf Kanäle zu montieren, ohne auf die Notwendigkeit von separaten Montageplatten.

Benötigen Sie Ihren Temperatursensor kalibriert?

Nach ISO17025 (UKAS) Kalibrierstandards rückverfolgbar, bietet Labfacility die Wahl zwischen einer 3- oder 5-Punkt-Rückverfolgbaren Kalibrierung mit unserer hauseigenen Kalibrieranlage, sie können beliebige Temperaturpunkte zwischen -10°C & +200°C auswählen.

Bitte beachten Sie, dass, wenn Sie eine Kalibrierungsoption zu Ihrer Bestellung hinzufügen, dies 2-3 Tage Vorlaufzeit zu Ihrer Bestellung hinzufügt.

Specifications

General Description	Compatible with Pt100 platinum resistance thermometers
Sensor Type	Surface Mounting for remote installation
Cable Length	1 or 2 metre
Cable Type	PFA screened or Fiberglass stainless steel overbraid available
Core / Strands	7/0.2mm
Head Type	High Impact ABS Plastic Construction
Thread	Entry Threads M16
Transmitter	Yes
Range	-50°C to +150°C or 0°C to +400°C (either can be re-ranged if desired)
Output	4-20 mA output