

**Handheld zwischen Packung / flacher
Lebensmittelsonde Typ T IEC 4,0 mm
Durchmesser x 140 mm Länge**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Runder Griff Stil Hand gehalten zwischen den Pack Sonden - Typ T IEC 4.0mm x 140mm

Specifications

Product Code	XE-3849-001
General Description	Hand Held Between Pack/Flat Food Probe
Thermocouple Type	T IEC
Thermocouple Junction	Single element, grounded junction offering fast response.
Probe Diameter	4.0mm
Probe Length	140mm approximately below the handle
Cable Length	2 metre (extended)
Cable Type	Polyurethane
Cable Termination	Moulded on miniature plug
Max. Temperature	+260°C
Min. Temperature	-50°C
Standards Met	IEC 584