

**Handheld zwischen Packung / flacher  
Lebensmittelsonde Typ T IEC 4,0 mm  
Durchmesser x 140 mm Länge**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Runder Griff Stil Hand gehalten zwischen den Pack Sonden - Typ T IEC 4.0mm x 140mm

## Specifications

<b>Product Code</b>	XE-3849-001
<b>General Description</b>	Hand Held Between Pack/Flat Food Probe
<b>Thermocouple Type</b>	T IEC
<b>Thermocouple Junction</b>	Single element, grounded junction offering fast response.
<b>Probe Diameter</b>	4.0mm
<b>Probe Length</b>	140mm approximately below the handle
<b>Cable Length</b>	2 metre (extended)
<b>Cable Type</b>	Polyurethane
<b>Cable Termination</b>	Moulded on miniature plug
<b>Max. Temperature</b>	+260°C
<b>Min. Temperature</b>	-50°C
<b>Standards Met</b>	IEC 584